


PHLDA3 Antibody (R35540)

Catalog No.	Formulation	Size
R35540-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23612
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PHLDA3 antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human PHLDA3-MYC probed with PHLDA3 antibody [1ug/ml] in Lane A and probed with anti-MYC [1/1000] in Lane C. Mock-transfected HEK293 probed with PHLDA3 [1ug/ml] in Lane B. Predicted molecular weight: ~14kDa.

Description

Additional name(s) for this target protein: Pleckstrin homology-like domain, family A, member 3; TIH1

Application Notes

Optimal dilution of the PHLDA3 antibody should be determined by the researcher.

Immunogen

Amino acids RARQSLGTGTLVS were used as the immunogen for this PHLDA3 antibody.

Storage

Aliquot and store the PHLDA3 antibody at -20oC.