


## Phospho-CD244 Antibody (R33769)

Catalog No.	Formulation	Size
R33769-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51744
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This phospho-CD244 antibody is available for research use only.



Western blot testing of immunoprecipitate from a lysate of human natural killer cells with (+) or without (-) pervanadate treatment of the cells using phospho-CD244 antibody at 2ug/ml.

## Description

Additional name(s) for this target protein: NAIL, NKR2B4, Nmrk, SLAMF4

## Application Notes

Optimal dilution of the phospho-CD244 antibody should be determined by the researcher.

## Immunogen

Amino acids EFLTIpYEDVKD were used as the immunogen for this phospho-CD244 antibody.

## **Storage**

Aliquot and store the phospho-CD244 antibody at -20oC.