


Phospho-Src Antibody (pY530) (F48453)


Catalog No.	Formulation	Size
F48453-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48453-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Chicken, Mouse, Rat, Zebrafish
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P12931
Applications	Dot Blot : 1:500 IHC (Paraffin) : 1:50-1:100
Limitations	This phospho-Src antibody is available for research use only.



Dot blot analysis of phospho-Src antibody. 50ng of phos-peptide or nonphos-peptide per dot were spotted.



IHC analysis of FFPE human hepatocarcinoma tissue stained with the phospho-Src antibody.

Description

This protein has tyrosine kinase activity and belongs to the src family of proteins.

Application Notes

Titration of the phospho-Src antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This phospho-Src antibody was produced from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding pY530 of human Src.

Storage

Aliquot the phospho-Src antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.