


# Phospho-TSC2 / Tuberin Antibody (pS664) (F48519)


Catalog No.	Formulation	Size
F48519-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48519-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P49815
Applications	Immunofluorescence : 1:10-1:50  Dot Blot : 1:500
Limitations	This phospho-Tuberin antibody is available for research use only.



Confocal immunofluorescent analysis of phospho-Tuberin antibody with MCF-7 cells followed by Alexa Fluor 488-conjugated goat anti-rabbit IgG (green). Actin filaments have been labeled with Alexa Fluor 555 Phalloidin (red).



Dot blot analysis of phospho-Tuberin antibody. 50ng of phos-peptide or nonphos-peptide per dot were spotted.

#### **Description**

Mutations in TSC2/Tuberin lead to tuberous sclerosis complex. The protein is believed to be a tumor suppressor and is able to specifically stimulate the intrinsic GTPase activity of the Ras-related protein RAP1A and RAB5. The protein associates with hamartin in a cytosolic complex, possibly acting as a chaperone for hamartin. TSC2 may have a function in vesicular transport, but may also play a role in the regulation of cell growth arrest and in the regulation of transcription mediated by steroid receptors. Interaction between TSC1 and TSC2 may facilitate vesicular docking.

## **Application Notes**

Titration of the phospho-Tuberin antibody may be required due to differences in protocols and secondary/substrate sensitivity.

#### **Immunogen**

This phospho-Tuberin antibody was produced from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding pS664 of human TSC2.

### **Storage**

Aliquot the phospho-Tuberin antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.