


## PHR1 Antibody (R35093)

Catalog No.	Formulation	Size
R35093-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	58473
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PHR1 antibody is available for research use only.



Western blot testing of rat brain lysate with PHR1 antibody at 1ug/ml. Predicted molecular weight: ~25/27kDa (rat/human).

## Description

Additional name(s) for this target protein: Pleckstrin homology domain containing family B member 1,; PLEKHB1

## Application Notes

Optimal dilution of the PHR1 antibody should be determined by the researcher.

1. This PHR1 antibody will detect isoforms a and c.

## Immunogen

Amino acids PPDSALESPFEE were used as the immunogen for this PHR1 antibody.

## Storage

Aliquot and store the PHR1 antibody at -20oC.