


PIAS2 Antibody (R34391)

Catalog No.	Formulation	Size
R34391-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9063
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PIAS2 antibody is available for research use only.



IHC testing of FFPE human testis with PIAS2 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Protein inhibitor of activated STAT, 2

Application Notes

Optimal dilution of the PIAS2 antibody should be determined by the researcher.

Immunogen

Amino acids KEAMKVSSQPCTKIE were used as the immunogen for this PIAS2 antibody.

Storage

Aliquot and store the PIAS2 antibody at -20oC.