


## PIAS2 Antibody (R35123)

Catalog No.	Formulation	Size
R35123-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9063
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PIAS2 antibody is available for research use only.



Western blot testing of human peripheral lymphocytes lysate with PIAS2 antibody at 1ug/ml. Predicted molecular weight: ~63kDa.

### Description

Additional name(s) for this target protein: Protein inhibitor of activated STAT 2

### Application Notes

Optimal dilution of the PIAS2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids PQQVREICISRD were used as the immunogen for this PIAS2 antibody.

## **Storage**

Aliquot and store the PIAS2 antibody at -20oC.