


PICK1 Antibody (R34330)

Catalog No.	Formulation	Size
R34330-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9463
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PICK1 antibody is available for research use only.



Western blot testing of mouse brain lysate with PICK1 antibody at 1ug/ml. Predicted molecular weight: ~47kDa.

Description

Additional name(s) for this target protein: Protein interacting with C kinase 1, Protein interacting with PRKCA 1

Application Notes

Optimal dilution of the PICK1 antibody should be determined by the researcher.

Immunogen

Amino acids HYNKLQADPKQ were used as the immunogen for this PICK1 antibody.

Storage

Aliquot and store the PICK1 antibody at -20oC.