


PIM2 Antibody (R35096)

Catalog No.	Formulation	Size
R35096-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11040
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PIM2 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human PIM2-MYC with PIM2 antibody (1ug/ml) in Lane A and anti-MYC (1/1000) in lane C. Mock-transfected HEK293 probed with PIM2 (1ug/ml) in Lane B. Predicted molecular weight: ~34kDa.

Description

Additional name(s) for this target protein: Serine/threonine protein kinase pim-2

Application Notes

Optimal dilution of the PIM2 antibody should be determined by the researcher.

Immunogen

Amino acids PPGGKDREAFEAE were used as the immunogen for this PIM2 antibody.

Storage

Aliquot and store the PIM2 antibody at -20oC.