


## PIM2 Antibody (R35096)

Catalog No.	Formulation	Size
R35096-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11040
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PIM2 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing human PIM2-MYC with PIM2 antibody (1ug/ml) in Lane A and anti-MYC (1/1000) in lane C. Mock-transfected HEK293 probed with PIM2 (1ug/ml) in Lane B. Predicted molecular weight: ~34kDa.

## Description

Additional name(s) for this target protein: Serine/threonine protein kinase pim-2

## Application Notes

Optimal dilution of the PIM2 antibody should be determined by the researcher.

## Immunogen

Amino acids PPGGKDREAFAEAE were used as the immunogen for this PIM2 antibody.

## **Storage**

Aliquot and store the PIM2 antibody at -20oC.