


Pink1 Antibody (F43671)

Catalog No.	Formulation	Size
F43671-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F43671-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q99MQ3
Applications	Western Blot : 1:1000
Limitations	This Pink1 antibody is available for research use only.



Description

Pink1 protects against mitochondrial dysfunction during cellular stress, potentially by phosphorylating mitochondrial proteins (By similarity).

Application Notes

Titration of the Pink1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 118-146 from the mouse protein was used as the immunogen for this Pink1 antibody.

Storage

Aliquot the Pink1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.