


## PINX Antibody (R35976)

Catalog No.	Formulation	Size
R35976-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54984
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PINX antibody is available for research use only.



Western blot testing of Jurkat lysate with PINX antibody at 2ug/ml. Predicted molecular weight: ~37kDa.

## Description

Additional name(s) for this target protein: PIN2-interacting protein 1, PINX1

## Application Notes

Optimal dilution of the PINX antibody should be determined by the researcher.

## Immunogen

Amino acids SMLAERRRKQKWAV were used as the immunogen for this PINX antibody.

## **Storage**

Aliquot and store the PINX antibody at -20oC.