


## Plakophilin 2 Antibody (R35623)

Catalog No.	Formulation	Size
R35623-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5318
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Plakophilin 2 antibody is available for research use only.



Western blot testing of mouse heart lysate with Plakophilin 2 antibody at 0.3ug/ml.  
Predicted molecular weight: 93-97 kDa.

### Description

Additional name(s) for this target protein: PKP2

### Application Notes

Optimal dilution of the Plakophilin 2 antibody should be determined by the researcher.

## **Immunogen**

Amino acids RGTAQYSSQKSVEER were used as the immunogen for this Plakophilin 2 antibody.

## **Storage**

Aliquot and store the Plakophilin 2 antibody at -20oC.