


Pleckstrin Antibody (R34045)

Catalog No.	Formulation	Size
R34045-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5341
Applications	Western Blot : 0.05-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Pleckstrin antibody is available for research use only.



Western blot testing of human PBMC lysate with Pleckstrin antibody at 0.05ug/ml.
Predicted molecular weight: ~40kDa, observed here at ~45kDa.

Description

Additional name(s) for this target protein: P47; PLEK

Application Notes

Optimal dilution of the Pleckstrin antibody should be determined by the researcher.

Immunogen

Amino acids EWIKAIQMASRTGK were used as the immunogen for this Pleckstrin antibody.

Storage

Aliquot and store the Pleckstrin antibody at -20oC.