


## PLEK Antibody / Pleckstrin [clone 28P96] (RQ8727)

Catalog No.	Formulation	Size
RQ8727	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

<b>Availability</b>	1-3 days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	28P96
<b>Purity</b>	Affinity chromatography
<b>UniProt</b>	P08567
<b>Applications</b>	Western Blot : 1:500-1:2000
<b>Limitations</b>	This PLEK antibody is available for research use only.



Western blot testing of human Raji cell lysate with PLEK antibody. Predicted molecular weight ~40 kDa.

### Description

Major protein kinase C substrate of platelets. [UniProt]

### Application Notes

Optimal dilution of the PLEK antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide specific to human Pleckstrin protein was used as the immunogen for the PLEK Antibody.

## **Storage**

Store the PLEK antibody at -20oC.