


## PLK3 Antibody (R34981)


Catalog No.	Formulation	Size
R34981-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1263
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PLK3 antibody is available for research use only.



Western blot testing of HeLa lysate with PLK3 antibody at 2ug/ml. Predicted molecular weight: ~71kDa.



IHC staining of FFPE human kidney with PLK3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Polo-like kinase 3

## Application Notes

Optimal dilution of the PLK3 antibody should be determined by the researcher.

## Immunogen

Amino acids PDRLPISSCVTPD were used as the immunogen for this PLK3 antibody.

## Storage

Aliquot and store the PLK3 antibody at -20oC.