


PLK3 Antibody (R34981)


Catalog No.	Formulation	Size
R34981-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1263
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PLK3 antibody is available for research use only.



Western blot testing of HeLa lysate with PLK3 antibody at 2ug/ml. Predicted molecular weight: ~71kDa.



IHC staining of FFPE human kidney with PLK3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Polo-like kinase 3

Application Notes

Optimal dilution of the PLK3 antibody should be determined by the researcher.

Immunogen

Amino acids PDRLPISSCVTPD were used as the immunogen for this PLK3 antibody.

Storage

Aliquot and store the PLK3 antibody at -20oC.