


PLK5 Antibody (R35955)

Catalog No.	Formulation	Size
R35955-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	126520
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PLK5 antibody is available for research use only.



Western blot testing of cerebral cortex lysate with PLK5 antibody at 2ug/ml.
Predicted/observed molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: Polo-like kinase 5

Application Notes

Optimal dilution of the PLK5 antibody should be determined by the researcher.

Immunogen

Amino acids SLDVPRSHGCAPTT were used as the immunogen for this PLK5 antibody.

Storage

Aliquot and store the PLK5 antibody at -20oC.