


PLRG1 Antibody (R35174)

Catalog No.	Formulation	Size
R35174-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5356
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PLRG1 antibody is available for research use only.



Western blot testing of A) human HEK293 and B) mouse NIH3T3 nuclear cell lysates with PLRG1 antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.

Description

Additional name(s) for this target protein: Pleiotropic regulator 1

Application Notes

Optimal dilution of the PLRG1 antibody should be determined by the researcher.

1. This PLRG1 antibody is specific for isoforms 1 and 2.

Immunogen

Amino acids PVSWKPEIHKRKF were used as the immunogen for this PLRG1 antibody.

Storage

Aliquot and store the PLRG1 antibody at -20oC.