


## PLUNC Antibody / SPLUNC1 / BPIFA1 (R34792)

Catalog No.	Formulation	Size
R34792-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51297
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PLUNC antibody is available for research use only.



Western blot testing of secretions from human primary airway cells in culture (lanes 1 and 2), and in human bronchoalveolar lavage fluid (lanes 3 and 4) with PLUNC antibody at 2ug/ml.

### Description

Additional name(s) for this target protein: short PLUNC1, SPURT; BPIFA1; BPI fold-containing family A member 1

### Application Notes

Optimal dilution of the PLUNC antibody should be determined by the researcher.

### Immunogen

Amino acids NEVLRGLDITLVHD were used as the immunogen for this PLUNC antibody.

## Storage

Aliquot and store the PLUNC antibody at -20oC.