


PLUNC Antibody / SPLUNC1 / BPIFA1 (R34792)

Catalog No.	Formulation	Size
R34792-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51297
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PLUNC antibody is available for research use only.



Western blot testing of secretions from human primary airway cells in culture (lanes 1 and 2), and in human bronchoalveolar lavage fluid (lanes 3 and 4) with PLUNC antibody at 2ug/ml.

Description

Additional name(s) for this target protein: short PLUNC1, SPURT; BPIFA1; BPI fold-containing family A member 1

Application Notes

Optimal dilution of the PLUNC antibody should be determined by the researcher.

Immunogen

Amino acids NEVLRGLDITLVHD were used as the immunogen for this PLUNC antibody.

Storage

Aliquot and store the PLUNC antibody at -20oC.