


PNP Antibody (R35905)

Catalog No.	Formulation	Size
R35905-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4860
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PNP antibody is available for research use only.



Western blot testing of human liver lysate with PNP antibody at 0.1ug/ml. Predicted molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: Purine nucleoside phosphorylase

Application Notes

Optimal dilution of the PNP antibody should be determined by the researcher.

Immunogen

Amino acids SGQNPLRGPND E were used as the immunogen for this PNP antibody.

Storage

Aliquot and store the PNP antibody at -20oC.