


Podocin Antibody (R36011)

Catalog No.	Formulation	Size
R36011-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Pig
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7827
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Podocin antibody is available for research use only.



kDa 200	Western blot of pig kidney lysate with Podocin antibody at 1ug/ml.
116 97	
66	
44	
31 –	
22 –	
14	
6-	
_	

Description

Additional name(s) for this target protein: PDCN, NPHS2, nephrosis 2

Application Notes

Optimal dilution of the Podocin antibody should be determined by the researcher.

Immunogen

Amino acids SPSKPVEPLNPKK were used as the immunogen for this Podocin antibody.

Storage

Aliquot and store the Podocin antibody at -20oC.