


Podocin Antibody (R36011)


Catalog No.	Formulation	Size
R36011-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Pig
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7827
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Podocin antibody is available for research use only.



Western blot testing of human [A] and mouse [B] kidney lysates with Podocin antibody at 0.3ug/ml. Predicted molecular weight ~42kDa.



Western blot of pig kidney lysate with Podocin antibody at 1ug/ml.

Description

Additional name(s) for this target protein: PDCN, NPHS2, nephrosis 2

Application Notes

Optimal dilution of the Podocin antibody should be determined by the researcher.

Immunogen

Amino acids SPSKPVEPLNPKK were used as the immunogen for this Podocin antibody.

Storage

Aliquot and store the Podocin antibody at -20oC.