


PODXL Antibody (R33559)

Catalog No.	Formulation	Size
R33559-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5420
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PODXL antibody is available for research use only.



Western blot testing of human kidney lysate with PODXL antibody at 1ug/ml. Predicted molecular weight: isoform 1 observed at ~59/165kDa (unmodified/glycosylated), isoform 2 observed at 55/160kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: Podocalyxin; Podocalyxin-like

Application Notes

Optimal dilution of the PODXL antibody should be determined by the researcher.

Immunogen

Amino acids DNLTKDDLDEEEDTH were used as the immunogen for this PODXL antibody.

Storage

Aliquot and store the PODXL antibody at -20oC.

References (1)