


POLDIP2 Antibody (R35608)

Catalog No.	Formulation	Size
R35608-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26073
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This POLDIP2 antibody is available for research use only.



Western blot testing of rat VSMC lysate with POLDIP2 antibody at 0.5ug/ml (Cat # R35068 and [R34668](#)). Signal is reduced by specific siRNA but not control.

Description

Additional name(s) for this target protein: Polymerase (DNA-directed), delta interacting protein 2

Application Notes

Optimal dilution of the POLDIP2 antibody should be determined by the researcher.

Immunogen

Amino acids RFERPDGSHFDVR were used as the immunogen for this POLDIP2 antibody.

Storage

Aliquot and store the POLDIP2 antibody at -20oC.