


## Poly (ADP-ribose) polymerase 1 Antibody / PARP1 [clone 1E12.] (RQ6272)


Catalog No.	Formulation	Size
RQ6272	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal (mouse origin)
<b>Isotype</b>	Mouse IgG2a
<b>Clone Name</b>	1E12.
<b>Purity</b>	Affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	P09874
<b>Localization</b>	Nuclear
<b>Applications</b>	Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 2-5ug/ml Immunofluorescence : 5ug/ml Flow Cytometry : 1-3ug/million cells
<b>Limitations</b>	This Poly (ADP-ribose) polymerase 1 antibody is available for research use only.




Poly (ADP-ribose) polymerase 1 Antibody Gallbladder Adenocarcinoma IHC.  
Immunohistochemistry staining of FFPE human gallbladder adenocarcinoma with Poly (ADP-ribose) polymerase 1 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.




Poly (ADP-ribose) polymerase 1 Antibody Ovarian Serous Adenocarcinoma IHC. Immunohistochemistry staining of FFPE human ovarian serous adenocarcinoma with Poly (ADP-ribose) polymerase 1 antibody. HIER: boil tissue sections in pH8 EDTA for 20 min and allow to cool before testing.



Poly (ADP-ribose) polymerase 1 Antibody WB. Western blot testing of human 1) MCF7, 2) HeLa, 3) Daudi, 4) HEK293, 5) rat testis and 6) mouse testis lysate with Poly (ADP-ribose) polymerase 1 antibody. Predicted molecular weight ~116 kDa.



Poly (ADP-ribose) polymerase 1 Antibody IF. Immunofluorescent staining of FFPE human A431 cells with Poly (ADP-ribose) polymerase 1 antibody (green) and DAPI nuclear stain (blue). HIER: steam section in pH6 citrate buffer for 20 min.



Poly (ADP-ribose) polymerase 1 Antibody FACS. Flow cytometry testing of human HL60 cells with Poly (ADP-ribose) polymerase 1 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= PARP antibody.

## Description

Poly [ADP-ribose] polymerase 1 Antibody detects PARP1, also known as ADPRT or PPOL, an enzyme that in humans is encoded by the PARP1 gene. PARP1 gene is mapped to 1q42.12. This gene encodes a chromatin-associated enzyme, poly(ADP-ribosyl)transferase, which modifies various nuclear proteins by poly(ADP-ribosyl)ation. The modification is dependent on DNA and is involved in the regulation of various important cellular processes such as differentiation, proliferation, and tumor transformation and also in the regulation of the molecular events involved in the recovery of cell from DNA damage. In addition, this enzyme may be the site of mutation in Fanconi anemia, and may participate in the pathophysiology of type I diabetes.

Additional PARP1 pathway and chromatin-associated DNA repair studies may benefit from our PARP1 antibody page featuring [recombinant rabbit monoclonal clone CFD-16](#) with knockdown-validated target recognition.

## Application Notes

Optimal dilution of the Poly (ADP-ribose) polymerase 1 antibody should be determined by the researcher.

## Immunogen

A human recombinant partial protein (amino acids Q670-R858) was used as the immunogen for the Poly (ADP-ribose) polymerase 1 antibody.

## Storage

After reconstitution, the Poly (ADP-ribose) polymerase 1 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.