


PON1 Antibody (R33755)

Catalog No.	Formulation	Size
R33755-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5444
Localization	Cytoplasmic, secreted
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PON1 antibody is available for research use only.



Western blot testing of human liver lysate with PON1 antibody at 1ug/ml. Predicted molecular weight: ~40 kDa.

Description

Additional name(s) for this target protein: Paraoxonase 1

Application Notes

Optimal dilution of the PON1 antibody should be determined by the researcher.

Immunogen

Amino acids EEPKVTQVYAEN were used as the immunogen for this PON1 antibody.

Storage

Aliquot and store the PON1 antibody at -20oC.