



POTE2 Antibody / POTEE (F54530)


| Catalog No. | Formulation | Size |
|---------------|--|---------|
| F54530-0.4ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.4 ml |
| F54530-0.08ML | In 1X PBS, pH 7.4, with 0.09% sodium azide | 0.08 ml |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|---|
| Species Reactivity | Human |
| Format | Purified |
| Clonality | Polyclonal (rabbit origin) |
| Isotype | Rabbit IgG |
| Purity | Antigen affinity purified |
| UniProt | Q6S8J3 |
| Applications | Western Blot : 1:500-1:2000 Flow Cytometry : 1:25 (1x10e6 cells) |
| Limitations | This POTE2 antibody is available for research use only. |







Flow cytometry testing of fixed and permeabilized human K562 cells with POTE2 antibody; Blue=isotype control, Green= POTE2 antibody.

Application Notes

The stated application concentrations are suggested starting points. Titration of the POTE2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 380-409 from the human protein was used as the immunogen for the POTE2 antibody.

Storage

Aliquot the POTE2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.