

POTE2 Antibody / POTE2 (F54530)


Catalog No.	Formulation	Size
F54530-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F54530-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Antigen affinity purified
UniProt	Q6S8J3
Applications	Western Blot : 1:500-1:2000 Flow Cytometry : 1:25 (1x10 ⁶ cells)
Limitations	This POTE2 antibody is available for research use only.



Western blot testing of human LNCaP cell lysate with POTE2 antibody. Predicted molecular weight ~121 kDa.



Western blot testing of human K562 cell lysate with POTE2 antibody. Predicted molecular weight ~121 kDa.



Flow cytometry testing of fixed and permeabilized human K562 cells with POTE2 antibody; Blue=isotype control, Green= POTE2 antibody.

Application Notes

The stated application concentrations are suggested starting points. Titration of the POTE2 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 380-409 from the human protein was used as the immunogen for the POTE2 antibody.

Storage

Aliquot the POTE2 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.