


## POU3F2 Antibody (R33161)

Catalog No.	Formulation	Size
R33161-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5454
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This POU3F2 antibody is available for research use only.



Western blot testing of mouse brain lysate with POU3F2 antibody at 0.1ug/ml. Predicted molecular weight: ~47kDa.

### Description

Additional name(s) for this target protein: POU class 3 homeobox 2

### Application Notes

Optimal dilution of the POU3F2 antibody should be determined by the researcher.

## Immunogen

Amino acids AQSLVQGDYGALQSN were used as the immunogen for this POU3F2 antibody.

## Storage

Aliquot and store the POU3F2 antibody at -20oC.