


Zebrafish / Danio rerio Pou3f3a Antibody (F52241)

Catalog No.	Formulation	Size
F52241-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52241-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P56224
Applications	Western Blot : 1:1000
Limitations	This Pou3f3a antibody is available for research use only.



Description

Transcription factor that may play important roles in patterning the embryonic brain.

This Zebrafish antibody is part of a [broader Zebrafish / Danio rerio antibody panel](#) offered by NSJ Bioreagents.

Application Notes

Titration of the Pou3f3a antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 381-407 from the zebrafish protein was used as the immunogen for this Pou3f3a antibody.

Storage

Aliquot the Pou3f3a antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.