


PPAR delta Antibody (R33249)

Catalog No.	Formulation	Size
R33249-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5467
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PPAR delta antibody is available for research use only.



Western blot testing of human MOLT4 lysate with PPAR delta antibody at 0.3ug/ml.
Predicted molecular weight: ~50 kDa.

Description

Additional name(s) for this target protein: Peroxisome proliferative activated receptor delta; PPARD

Application Notes

Optimal dilution of the PPAR delta antibody should be determined by the researcher.

1. This PPAR delta antibody is specific for isoform 1 only.

Immunogen

Amino acids CHPLLQEIYKDMY were used as the immunogen for this PPAR delta antibody.

Storage

Aliquot and store the PPAR delta antibody at -20oC.