


## PPP3CC Antibody (F40087)


Catalog No.	Formulation	Size
F40087-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40087-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Purified
<b>UniProt</b>	P48454
<b>Applications</b>	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
<b>Limitations</b>	This PPP3CC antibody is available for research use only.



Western blot analysis of PPP3CC antibody and HL-60, Jurkat, A375 lysate.



IHC analysis of FFPE human breast carcinoma stained with the PPP3CC antibody

## Description

Calmodulin-dependent protein phosphatase, calcineurin, is involved in a wide range of biologic activities, acting as a Ca(2+)-dependent modifier of phosphorylation status. In testis, the motility of the sperm is thought to be controlled by cAMP-dependent phosphorylation and a unique form of calcineurin appears to be associated with the flagellum. The calcineurin holoenzyme is composed of catalytic and regulatory subunits of 60 and 18 kD, respectively. At least 3 genes, calcineurin A-alpha, calcineurin A-beta, and calcineurin A-gamma (CALNA3), have been cloned for the catalytic subunit. These genes have been identified in humans, mice, and rats, and are highly conserved between species (90 to 95% amino acid identity).

## Application Notes

Titration of the PPP3CC antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 18-47 from the human protein was used as the immunogen for this PPP3CC antibody.

## Storage

Aliquot the PPP3CC antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.