

PRAM1 Antibody (F53937)

Catalog No.	Formulation	Size
F53937-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q96QH2
Applications	Western Blot : 1:500-1000
Limitations	This PRAM1 antibody is available for research use only.



Western blot testing of human THP-1 cell lysate with PRAM1 antibody at 1:1000.
Predicted molecular weight: 79 kDa but can also be observed at ~97 kDa.

Description

PML-RARA-regulated adapter molecule 1 may be involved in myeloid differentiation. May be involved in integrin signaling in neutrophils. Binds to PtdIns(4)P.

Application Notes

Titration of the PRAM1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 587-618 from the human protein was used as the immunogen for the PRAM1 antibody.

Storage

Aliquot the PRAM1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.