


## PRAM1 Antibody (R33436)

Catalog No.	Formulation	Size
R33436-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84106
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PRAM1 antibody is available for research use only.



Western blot testing of Jurkat lysate with PRAM1 antibody at 2ug/ml. Predicted molecular weight: ~79kDa but routinely observed at 97kDa.

## Description

Additional name(s) for this target protein: PML-RARA target gene encoding an Adaptor Molecule-1

## Application Notes

Optimal dilution of the PRAM1 antibody should be determined by the researcher.

## **Immunogen**

Amino acids DFCDPLENQPLPLGR were used as the immunogen for this PRAM1 antibody.

## **Storage**

Aliquot and store the PRAM1 antibody at -20oC.