


PRDM14 Antibody (R33633)

Catalog No.	Formulation	Size
R33633-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	63978
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This PRDM14 antibody is available for research use only.



Western blot testing of A549 lysate with PRDM14 antibody at 2ug/ml. Predicted molecular weight: ~64kDa.

Description

Additional name(s) for this target protein: PR domain containing 14

Application Notes

Optimal dilution of the PRDM14 antibody should be determined by the researcher.

Immunogen

Amino acids DSHVRRSHKEDD were used as the immunogen for this PRDM14 antibody.

Storage

Aliquot and store the PRDM14 antibody at -20oC.