


PRDM4 Antibody / PR domain zinc finger protein 4 (R33335)

Catalog No.	Formulation	Size
R33335-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	11108
Applications	Western Blot : 0.3-2ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PRDM4 antibody is available for research use only.



Western blot testing of human brain lysate with PRDM4 antibody at 0.3ug/ml. Predicted molecular weight: ~88 kDa, observed here at ~70 kDa.

Description

Additional name(s) for this target protein: PR domain containing 4; PFM1

Application Notes

Optimal dilution of the PRDM4 antibody should be determined by the researcher.

Immunogen

Amino acids CSAVYSADESLSAHK were used as the immunogen for this PRDM4 antibody.

Storage

Aliquot and store the PRDM4 antibody at -20oC.