


## Prdm9 Antibody (R35239)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35239-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| Availability       | 1-3 business days                                         |
| Species Reactivity | Rat                                                       |
| Format             | Antigen affinity purified                                 |
| Host               | Goat                                                      |
| Clonality          | Polyclonal (goat origin)                                  |
| Isotype            | Goat Ig                                                   |
| Purity             | Antigen affinity                                          |
| Gene ID            | 365155                                                    |
| Applications       | Western Blot : 2-3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This Prdm9 antibody is available for research use only.   |



Western blot testing of rat ovary lysate with Prdm9 antibody at 2ug/ml. Predicted molecular weight: ~93kDa.

## Description

Additional name(s) for this target protein: PR domain containing 9

## Application Notes

Optimal dilution of the Prdm9 antibody should be determined by the researcher.

## Immunogen

Amino acids QEHFDSQNKNDK were used as the immunogen for this Prdm9 antibody.

## **Storage**

Aliquot and store the Prdm9 antibody at -20oC.