


## Prefoldin Antibody (R33557)

Catalog No.	Formulation	Size
R33557-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5201
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Prefoldin antibody is available for research use only.



Western blot testing of K562 lysate with Prefoldin antibody at 1ug/ml. Predicted molecular weight: ~14kDa. Both observed bands are blocked by addition of immunizing peptide.

## Description

Additional name(s) for this target protein: PFDN1

## Application Notes

Optimal dilution of the Prefoldin antibody should be determined by the researcher.

## **Immunogen**

Amino acids DNIREMLMARRAQ were used as the immunogen for this Prefoldin antibody.

## **Storage**

Aliquot and store the Prefoldin antibody at -20oC.