


PREX1 Antibody (R33979)


| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33979-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat, Dog, Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 57580 |
| Localization | Cytoplasmic |
| Applications | Western Blot : 1-3ug/ml IHC (FFPE) : 10ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This PREX1 antibody is available for research use only. |



Western blot testing of human temporal cortex lysate with PREX1 antibody at 1ug/ml.
Predicted molecular weight: ~186 kDa, observed here at ~150 kDa.



IHC testing of FFPE human testis tissue with PREX1 antibody at 10ug/ml. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

Description

Additional name(s) for this target protein: Phosphatidylinositol-3,4,5-trisphosphate-dependent Rac exchange factor 1

Application Notes

Optimal dilution of the PREX1 antibody should be determined by the researcher.

Immunogen

Amino acids ERSNLPTDASTTA were used as the immunogen for this PREX1 antibody.

Storage

Aliquot and store the PREX1 antibody at -20oC.