


PRMT2 Antibody (R33600)


Catalog No.	Formulation	Size
R33600-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3275
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This PRMT2 antibody is available for research use only.



Western blot testing of HeLa lysate with PRMT2 antibody at 2ug/ml. Predicted molecular weight: ~37kDa.



Western blot of mouse (A) and rat (B) brain lysates with PRMT2 antibody at 2ug/ml. Predicted molecular weight: ~37kDa.



Western blot of pig brain lysate with PRMT2 antibody at 2ug/ml. Predicted molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: Protein arginine methyltransferase 2

Application Notes

Optimal dilution of the PRMT2 antibody should be determined by the researcher.

Immunogen

Amino acids DQPRTTYHHSVILQ were used as the immunogen for this PRMT2 antibody.

Storage

Aliquot and store the PRMT2 antibody at -20oC.