


## PRMT3 Antibody (R33574)

Catalog No.	Formulation	Size
R33574-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	10196
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This PRMT3 antibody is available for research use only.



Western blot testing of MCF7 lysate with PRMT3 antibody at 1ug/ml. Predicted molecular weight: ~60kDa, routinely observed at 60~70kDa.

### Description

Additional name(s) for this target protein: Protein arginine methyltransferase 3

### Application Notes

Optimal dilution of the PRMT3 antibody should be determined by the researcher.

### Immunogen

Amino acids DPKTLISEPCGIKH were used as the immunogen for this PRMT3 antibody.

## **Storage**

Aliquot and store the PRMT3 antibody at -20oC.