


## PRMT5 Antibody (R35528)


Catalog No.	Formulation	Size
R35528-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10419
Applications	Immunofluorescence : 10ug/ml Flow Cytometry : 10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PRMT5 antibody is available for research use only.



IF/ICC staining of fixed and permeabilized human U-2 OS cells with PRMT5 antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



IF/ICC staining of fixed and permeabilized human A431 cells with PRMT5 antibody (green) at 10ug/ml and DAPI nuclear stain (blue).



FACS testing of fixed and permeabilized human A431 cells with PRMT5 antibody (blue) at 10ug/ml and naive goat Ig (black).

## Description

Additional name(s) for this target protein: Protein arginine methyltransferase 5

## Application Notes

Optimal dilution of the PRMT5 antibody should be determined by the researcher.

## Immunogen

Amino acids RACREKDRDPEAQ were used as the immunogen for this PRMT5 antibody.

## Storage

Aliquot and store the PRMT5 antibody at -20oC.