


PRMT6 Antibody (R33955)

Catalog No.	Formulation	Size
R33955-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55170
Applications	Western Blot : 0.5-2ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PRMT6 antibody is available for research use only.



Western blot testing of MCF7 lysate with PRMT6 antibody at 0.5ug/ml. Predicted molecular weight: ~42kDa.

Description

Additional name(s) for this target protein: Protein arginine methyltransferase 6

Application Notes

Optimal dilution of the PRMT6 antibody should be determined by the researcher.

Immunogen

Amino acids ERDQLYYECYSD were used as the immunogen for this PRMT6 antibody. This sequence is specific to isoform 1.

Storage

Aliquot and store the PRMT6 antibody at -20oC.