


## PRMT7 Antibody (R35135)

Catalog No.	Formulation	Size
R35135-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54496
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PRMT7 antibody is available for research use only.



Western blot testing of HeLa lysate with PRMT7 antibody at 0.3ug/ml. Predicted molecular weight: ~78kDa.

### Description

Additional name(s) for this target protein: Protein arginine methyltransferase 7

### Application Notes

Optimal dilution of the PRMT7 antibody should be determined by the researcher.

## **Immunogen**

Amino acids PRFGEINDQDRTDR were used as the immunogen for this PRMT7 antibody.

## **Storage**

Aliquot and store the PRMT7 antibody at -20oC.