


PRODH Antibody (R33952)

Catalog No.	Formulation	Size
R33952-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5625
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PRODH antibody is available for research use only.



IHC testing of FFPE human cerebellum with PRODH antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Proline dehydrogenase

Application Notes

Optimal dilution of the PRODH antibody should be determined by the researcher.

Immunogen

Amino acids ERDGS GTNKRDKQYQ were used as the immunogen for this PRODH antibody.

Storage

Aliquot and store the PRODH antibody at -20oC.