


## Progranulin Antibody (R35108)

Catalog No.	Formulation	Size
R35108-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2896
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Progranulin antibody is available for research use only.



Western blot testing of A431 lysate with Progranulin antibody at 1ug/ml. Predicted molecular weight: 64/70~88kDa (unmodified/glycosylated), observed here at ~75kDa.

## Description

Additional name(s) for this target protein: Granulin; GRN

## Application Notes

Optimal dilution of the Progranulin antibody should be determined by the researcher.

## Immunogen

Amino acids PQDTVCDLIQSK were used as the immunogen for this Progranulin antibody.

## Storage

Aliquot and store the Progranulin antibody at -20oC.