


Progranulin Antibody (R35108)

Catalog No.	Formulation	Size
R35108-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2896
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Progranulin antibody is available for research use only.



Western blot testing of A431 lysate with Progranulin antibody at 1ug/ml. Predicted molecular weight: 64/70~88kDa (unmodified/glycosylated), observed here at ~75kDa.

Description

Additional name(s) for this target protein: Granulin; GRN

Application Notes

Optimal dilution of the Progranulin antibody should be determined by the researcher.

Immunogen

Amino acids PQDTVCDLIQSK were used as the immunogen for this Progranulin antibody.

Storage

Aliquot and store the Progranulin antibody at -20oC.