


PrP Antibody (R35863)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35863-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Pig |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 5621 |
| Localization | Cytoplasmic |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This PrP antibody is available for research use only. |



Western blot testing of human brain lysate with PrP antibody at 0.3ug/ml. Expected molecular weight: 20~29kDa.

Description

Additional name(s) for this target protein: Prion protein, PRNP

Application Notes

Optimal dilution of the PrP antibody should be determined by the researcher.

Immunogen

Amino acids SDYEDRYREN were used as the immunogen for this PrP antibody.

Storage

Aliquot and store the PrP antibody at -20oC.